


X-Ray Fluorescence Spectrum Viewer

Process Data Plot Spectrum  Plot

Choose Spectra


Browse... No file selected

Project Name

Element:

(Fe) Iron ▼

First, load data, must be v24 CSV format




X-Ray Fluorescence Spectrum Viewer

Process Data Plot Spectrum Plot

Choose Spectra
Browse... No file selected

Project Name

Element:
(Fe) Iron



There is no limit, here we will use ~2,000 files

X-Ray Fluorescence Spectrum Viewer

Process Data Plot Spectrum Plot


Choose Spectra

Browse... 2026 files

Upload complete

Project Name


Element: (Fe) Iron



You should see a little window show up down here

Processing Data

X-Ray Fluorescence Spectrum Viewer

Process Data Plot Spectrum  Plot

Choose Spectra

Browse... 2026 files

Upload complete


Project Name

RSP15B

Element:

(Fe) Iron

You can add a name to the project here to help streamline the saving of plots



Enter Values [Table](#)

Elemental lines to show:

- Spectrum
- Ne.K.alpha
- Na.K.alpha
- Mg.K.alpha
- Al.K.alpha
- Si.K.alpha
- P.K.alpha
- S.K.alpha
- Cl.K.alpha
- Ar.K.alpha
- K.K.alpha
- Ca.K.alpha
- Sc.K.alpha
- Ti.K.alpha
- V.K.alpha
- Cr.K.alpha
- Mn.K.alpha
- Fe.K.alpha
- Co.K.alpha
- Ni.K.alpha
- Cu.K.alpha
- Zn.K.alpha
- Ga.K.alpha
- Ge.K.alpha
- As.K.alpha
- Se.K.alpha
- Br.K.alpha
- Kr.K.alpha
- Rb.K.alpha
- Sr.K.alpha
- Y.K.alpha
- Zr.K.alpha
- Nb.K.alpha
- Mo.K.alpha
- Mo.L.alpha
- Mo.L.beta

After some time, you will see a table show up in the Counts tab


Spectral Lines [Add Categories](#)

Show entries

Search:

	Spectrum	Ca.K.alpha	Ti.K.alpha	Fe.K.alpha	Cu.K.alpha	Zn.K.alpha	Pb.L.alpha
1	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00309.3mm.CSV	959.693183052729	40.7689085390009	428.098268071482	142.237825667004	140.968433852405	318.370061344596
2	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00311.3mm.CSV	897.427041705456	40.0105705718379	427.867469559737	129.840648464687	133.038856413157	303.87921192859
3	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00313.3mm.CSV	480.176303686001	45.2694795180334	443.314484238687	133.632338300502	137.951567591735	318.86462958405
4	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00315.3mm.CSV	510.625221628392	46.4894145086868	442.589117487488	132.874000333339	134.110420931974	313.902461581527
5	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00317.3mm.CSV	497.799418618549	43.1593216963625	408.233110453409	132.016748718285	136.270035577591	312.847382670692
6	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00319.3mm.CSV	440.000877034345	41.4942752902003	414.695468782276	133.79719438032	136.566776521263	307.423617644678
7	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00321.3mm.CSV	399.462767007089	45.7310765415239	465.833824741831	133.764223164356	135.511697610428	309.171092090749
8	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00323.3mm.CSV	414.893296078057	46.1267311330872	526.368977251015	134.423647483629	135.973294633918	310.81965288893
9	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00325.3mm.CSV	424.669261611267	50.0832770487201	578.628354553333	136.962431112826	136.319492401536	314.759713196581
10	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00327.3mm.CSV	421.421596838852	49.2095398256845	588.338377654616	135.000643762992	137.242686448517	321.485841253157

Showing 1 to 10 of 2,026 entries

Enter Values  Table


Elemental lines to

- Spectrum
- Ne.K.alpha
- Na.K.alpha
- Mg.K.alpha
- Al.K.alpha
- Si.K.alpha
- P.K.alpha
- S.K.alpha
- Cl.K.alpha
- Ar.K.alpha
- K.K.alpha
- Ca.K.alpha
- Sc.K.alpha
- Ti.K.alpha
- V.K.alpha
- Cr.K.alpha
- Mn.K.alpha
- Fe.K.alpha
- Co.K.alpha
- Ni.K.alpha
- Cu.K.alpha
- Zn.K.alpha
- Ga.K.alpha
- Ge.K.alpha
- As.K.alpha
- Se.K.alpha
- Br.K.alpha
- Kr.K.alpha
- Rb.K.alpha
- Sr.K.alpha
- Y.K.alpha
- Zr.K.alpha
- Nb.K.alpha
- Mo.K.alpha
- Mo.L.alpha
- Mo.L.beta

Go to "Add Categories" to augment the data with other information

Spectral Lines Add Categories

	Spectrum	Qualitative	Depth
1	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00309.3mm.CSV	a	0
	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00311.3mm.CSV	b	0
3	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00313.3mm.CSV	c	0
4	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00315.3mm.CSV	NULL	0
5	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00317.3mm.CSV	NULL	0
6	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00319.3mm.CSV	NULL	0
7	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00321.3mm.CSV	NULL	0
8	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00323.3mm.CSV	NULL	0
9	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00325.3mm.CSV	NULL	0
10	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00327.3mm.CSV	NULL	0
11	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00329.3mm.CSV	NULL	0
12	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00331.3mm.CSV	NULL	0
13	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00333.3mm.CSV	NULL	0
14	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00335.3mm.CSV	NULL	0
15	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00337.3mm.CSV	NULL	0
16	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00339.3mm.CSV	NULL	0
17	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00341.3mm.CSV	NULL	0
18	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00343.3mm.CSV	NULL	0
19	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00345.3mm.CSV	NULL	0
20	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00347.3mm.CSV	NULL	0
21	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00349.3mm.CSV	NULL	0
22	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00351.3mm.CSV	NULL	0
23	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00353.3mm.CSV	NULL	0
24	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00355.3mm.CSV	NULL	0
25	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00357.3mm.CSV	NULL	0
26	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00359.3mm.CSV	NULL	0
27	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00361.3mm.CSV	NULL	0
28	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00363.3mm.CSV	NULL	0
29	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00365.3mm.CSV	NULL	0
30	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00367.3mm.CSV	NULL	0
31	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00369.3mm.CSV	NULL	0
32	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00371.3mm.CSV	NULL	0
33	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00373.3mm.CSV	NULL	0
34	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00375.3mm.CSV	NULL	0
35	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00377.3mm.CSV	NULL	0
36	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH002CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305374.5mm.CSV	NULL	0
37	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH004CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305372.5mm.CSV	NULL	0
38	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH006CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305370.5mm.CSV	NULL	0
39	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH008CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305368.5mm.CSV	NULL	0
40	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH010CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305366.5mm.CSV	NULL	0
41	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH012CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305364.5mm.CSV	NULL	0

Enter Values  Table

Elemental lines to show:

- Spectrum
- Ne.K.alpha
- Na.K.alpha
- Mg.K.alpha
- Al.K.alpha
- Si.K.alpha
- P.K.alpha
- S.K.alpha
- Cl.K.alpha
- Ar.K.alpha
- K.K.alpha
- Ca.K.alpha
- Sc.K.alpha
- Ti.K.alpha
- V.K.alpha
- Cr.K.alpha
- Mn.K.alpha
- Fe.K.alpha
- Co.K.alpha
- Ni.K.alpha
- Cu.K.alpha
- Zn.K.alpha
- Ga.K.alpha
- Ge.K.alpha
- As.K.alpha
- Se.K.alpha
- Br.K.alpha
- Kr.K.alpha
- Rb.K.alpha
- Sr.K.alpha
- Y.K.alpha
- Zr.K.alpha
- Nb.K.alpha
- Mo.K.alpha
- Mo.L.alpha
- Mo.L.beta

Spectral Lines

	Spectrum	Qualitative	Depth
1	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00309.3mm.CSV	a	6
2	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00311.3mm.CSV	b	8
3	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00313.3mm.CSV	c	10
4	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00315.3mm.CSV	NULL	12
5	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00317.3mm.CSV	NULL	14
6	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00319.3mm.CSV	NULL	16
7	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00321.3mm.CSV	NULL	18
8	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00323.3mm.CSV	NULL	20
9	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00325.3mm.CSV	NULL	22
10	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00327.3mm.CSV	NULL	24
11	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00329.3mm.CSV	NULL	26
12	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00331.3mm.CSV	NULL	28
13	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00333.3mm.CSV	NULL	30
14	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00335.3mm.CSV	NULL	32
15	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00337.3mm.CSV	NULL	34
16	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00339.3mm.CSV	NULL	36
17	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00341.3mm.CSV	NULL	38
18	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00343.3mm.CSV	NULL	40
19	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00345.3mm.CSV	NULL	42
20	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00347.3mm.CSV	NULL	44
21	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00349.3mm.CSV	NULL	46
22	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00351.3mm.CSV	NULL	48
23	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00353.3mm.CSV	NULL	50
24	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00355.3mm.CSV	NULL	52
25	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00357.3mm.CSV	NULL	54
26	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00359.3mm.CSV	NULL	56
27	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00361.3mm.CSV	NULL	58
28	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00363.3mm.CSV	NULL	60
29	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00365.3mm.CSV	NULL	62
30	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00367.3mm.CSV	NULL	64
31	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00369.3mm.CSV	NULL	66
32	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00371.3mm.CSV	NULL	68
33	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00373.3mm.CSV	NULL	70
34	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00375.3mm.CSV	NULL	72
35	CYCLE_C1 COREID_RSP15B He 0-0073 DEPTH_00377.3mm.CSV	NULL	74
36	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH002CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305374.5mm.CSV	NULL	426
37	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH004CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305372.5mm.CSV	NULL	428
38	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH006CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305370.5mm.CSV	NULL	430
39	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH008CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305368.5mm.CSV	NULL	432
40	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH010CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305366.5mm.CSV	NULL	434
41	CYCLE_C1 COREID_RSP15B He 1000-0120 DEPTH012CYCLE_C1 COREID_RSP15-B He 1000-0120 DEPTH_305364.5mm.CSV	NULL	436

Here, I've added the depths for each assay

Enter Dates Run Age Model Age Results Table

Age Input

Age Curve

Age Table

	Depth	14C Age	Sigma	CalCurve
1	0.00	0.00	0.00	intcal13
2	0.00	0.00	0.00	intcal13
3	0.00	0.00	0.00	intcal13
4	0.00	0.00	0.00	intcal13
5	0.00	0.00	0.00	intcal13
6	0.00	0.00	0.00	intcal13
7	0.00	0.00	0.00	intcal13
8	0.00	0.00	0.00	intcal13
9	0.00	0.00	0.00	intcal13
10	0.00	0.00	0.00	intcal13
11	0.00	0.00	0.00	intcal13
12	0.00	0.00	0.00	intcal13
13	0.00	0.00	0.00	intcal13
14	0.00	0.00	0.00	intcal13
15	0.00	0.00	0.00	intcal13
16	0.00	0.00	0.00	intcal13
17	0.00	0.00	0.00	intcal13
18	0.00	0.00	0.00	intcal13
19	0.00	0.00	0.00	intcal13
20	0.00	0.00	0.00	intcal13

Next, go to the Age Model tab

Enter Dates Run Age Model Age Results Table


Age Input Age Curve Age Table

	Depth	14C Age	Sigma	CalCurve
1	405.00	1698.00	23.00	intcal13
2	605.00	4024.00	24.00	intcal13
3	805.00	5804.00	28.00	intcal13
4	1135.00	7942.00	27.00	intcal13
5	1405.00	6754.00	30.00	intcal13
6	2555.00	19364.00	47.00	intcal13
7	2975.00	19268.00	49.00	intcal13
8	4655.00	23251.00	73.00	intcal13
9	4655.00	23654.00	62.00	intcal13
10	4715.00	22037.00	94.00	intcal13
11	4715.00	23017.00	71.00	intcal13
12	0.00	0.00	0.00	intcal13
13	0.00	0.00	0.00	intcal13
14	0.00	0.00	0.00	intcal13
15	0.00	0.00	0.00	intcal13
16	0.00	0.00	0.00	intcal13
17	0.00	0.00	0.00	intcal13
18	0.00	0.00	0.00	intcal13
19	0.00	0.00	0.00	intcal13
20	0.00	0.00	0.00	intcal13

Add the depth, 14C age, and sigma for each date. You can choose a calibration curve to the right

Enter Dates Run Age Model Age Results Table

Age Input Age Curve Age Table



After some time, the Age Curve will appear in the named tab

K-Means

3

Colour

Cluster


Point Size

2

Elipse

[Plot](#) [Results](#)

Go to PCA next to perform principle components analysis and/or k-means cluster analysis



K-Means

3

Colour

Climate

Black

Cluster


Climate

Qualitative

Quantitative

Plot Results

You can color the points based on other values, like climatic period, cluster analysis, or other data you have added



K-Means

3

Colour

Climate


Point Size

2 15

Ellipse

Plot Results

You can also adjust the size of the points for clarity. Smaller for big data sets, bigger for small data sets



K-Means

3

Colour

Climate


Point Size

2 15

Elipse

[Plot](#) [Results](#)

Download the plot, and you will get a publication-ready tiff with the project name + plot type



K-Means

3

Colour

Climate


Point Size

2 15

Ellipse

Plot Results

You can add ellipses to the data as well



- Climatic Period**
- 1. Late Holocene
 - 2. Altithermal
 - 3. Holocene Transition
 - 4. Younger Dryas
 - 5. Bølling-Allerød
 - 6. Deglaciation
 - 7. Last Glacial Maximum
 - 8. Glacial

Create Plot

1 2 3 4 5

↓ 1 ↓ 2 ↓ 3 ↓ 4 ↓ 5

Element:
Fe.K.alpha

Ratio:
None

Time Series Type
Black

Smoothed Mean Average
1 50


Line Size
1 15

Point Size
1 5 15

X Axis
Depth

mm
 cm
 m
 inches

On the Timeseries tab, you can plot the data along the core's length



Create Plot

1 2 3 4 5

↓ 1 ↓ 2 ↓ 3 ↓ 4 ↓ 5

Element:
Fe.K.alpha

Ratio:
None

Time Series Type
Climate


- Black
- Smooth
- Ramp
- Cluster
- Climate
- Qualitative Point
- Qualitative Line
- Quantitative

Point Size
1 5 15

X Axis
Depth

- mm
- cm
- m
- inches

RSP15B_PCAPlot.tiff



As before, you can color based on Climate, Cluster Analysis, or a Qualitative attribute

Create Plot

1 2 3 4 5

↓ 1 ↓ 2 ↓ 3 ↓ 4 ↓ 5

Element:
Fe.K.alpha

Ratio:
Ca.K.alpha

- K.K.alpha
- Ca.K.alpha
- Sc.K.alpha
- Ti.K.alpha
- V.K.alpha
- Cr.K.alpha
- Mn.K.alpha
- Fe.K.alpha


Line Size
1 15

Point Size
1 5

X Axis
Depth

- mm
- cm
- m
- inches

RSP15B_PCAPlot.tiff



You can change elements, or even ratio them together by selecting them from the dropdown list (in order of atomic number)

Create Plot

1 2 3 4 5

↓ 1 ↓ 2 ↓ 3 ↓ 4 ↓ 5

Element:
Fe.K.alpha

Ratio:
Ca.K.alpha

Time Series Type
Climate

Smoothed Mean Average
1 50


Line Size
1 15

Point Size
1 5 15

X Axis
Age
Depth
Age

m
 inches

RSP15B_PCAPlot.tiff



Change the x-axis from Depth to Age (using the age model) here

Create Plot

1 2 3 4 5

↓ 1 ↓ 2 ↓ 3 ↓ 4 ↓ 5

Element:
Fe.K.alpha

Ratio:
Ca.K.alpha

Time Series Type
Climate

Smoothed Mean Average
1 10 50


Line Size
1 15

Point Size
1 5 15

X Axis
Age

mm
 cm
 m
 inches

Add a smoothed-mean average here



- Climatic Period
- 1. Late Holocene
 - 2. Altithermal
 - 3. Holocene Transition
 - 4. Younger Dryas
 - 5. Bølling-Allerød
 - 6. Deglaciation
 - 7. Last Glacial Maximum
 - 8. Glacial

Create Plot

1 2 3 4 5

↓ 1 ↓ 2 ↓ 3 ↓ 4 ↓ 5

Element:
Mg.K.alpha

Ratio:
Ca.K.alpha

Time Series Type
Climate

Smoothed Mean Average
1 3 50


Line Size
1 15

Point Size
1 5 15

X Axis
Age

mm
 cm
 m
 inches

You can work on up to 5 plots simultaneously



Ternary Diagrams are also available

Colour
Cluster

Axis A
Al.K.alpha


Axis B
Si.K.alpha

Axis C
Ca.K.alpha

Density Contour
 Normalize

Point Size
2 5 15

Plot



Cluster
1
2
3

Colour

Climate

- Black
- Cluster
- Climate
- Qualitative
- Quantitative

Si.K.alpha

Axis C

Ca.K.alpha

Density Contour


Normalize

Point Size

2 15

Plot

As before, you can color points by different data types



- Climatic Period
- 1. Late Holocene
 - 2. Altithermal
 - 3. Holocene Transition
 - 4. Younger Dryas
 - 5. Bølling-Allerød
 - 6. Deglaciation
 - 7. Last Glacial Maximum
 - 8. Glacial

Colour: Climate


Axis A: Na.K.alpha

Axis B: K.K.alpha

- Si.K.alpha
- P.K.alpha
- S.K.alpha
- Cl.K.alpha
- Ar.K.alpha
- K.K.alpha
- Ca.K.alpha
- Sr.K.alpha

Point Size: 2 (range 2 to 15)

Plot



- Climatic Period
- 1. Late Holocene
 - 2. Altithermal
 - 3. Holocene Transition
 - 4. Younger Dryas
 - 5. Bølling-Allerød
 - 6. Deglaciation
 - 7. Last Glacial Maximum
 - 8. Glacial

Each axis can be changed to your preference as well

Ratio Plot Type
Climate

Element A
Na.K.alpha

Element B
K.K.alpha


Element C
Mg.K.alpha

Element D
Ca.K.alpha

Point Size
2 3 15

Elipse

Plot



- Climatic Period
- 1. Late Holocene
 - 2. Altithermal
 - 3. Holocene Transition
 - 4. Younger Dryas
 - 5. Bølling-Allerød
 - 6. Deglaciation
 - 7. Last Glacial Maximum
 - 8. Glacial