

# Exploring Russian-American Trade through Comparison of Chemical XRF Signatures of Glass from Colonial Russian Sites in Alaska and the Tal'tsinka Factory in Central Siberia

Dr. Claudia Brackett, California State University Stanislaus  
Julia Kleyman, Thermo Fisher Scientific  
David McMahan, Alaska Office of History Archaeology










Image NASA  
Image © 2008 TerraMetrics

©2007 Google™

# Иркутск








1725: Bering's trek across Siberia from Irkutsk to Okhotsk (in red)

# Tal'tsinka Glassworks: 1784

Making Glass III




(Generalized illustration from Diderot)

# Archaeology at the Baranov-Laxman Glass Factory: 2004-









SIT-00002

SIT-00002

**NOOW TLEIN - CASTLE HILL**  
ARCHAEOLOGICAL PROJECT




**Sitka, Alaska 1998**

ADNR Office of History & Archaeology

# Castle Hill Workshop Area / Architecture








Von Kittlitz 1827





LOCATION	INSPECTOR	COR 1	NOTE	Mo	Zr	Sr	U
Alaska-Inkutzk	CLB	IKT-04-1		0.81	42.096	82.255	1.24
Alaska-Inkutzk	CLB	IKT-04-1		0.27	Not available	9.612	0.75
Alaska-Inkutzk	CLB	IKT-04-2		0.12	114.151	122.142	26.285
Alaska-Inkutzk	CLB	IKT-04-2		2.366	1.016	Not available	1.86
Alaska-Inkutzk	CLB	IKT-04-3		0.5	95.935	105.319	1.09
Alaska-Inkutzk	CLB	IKT-04-3		0.45	Not available	14.814	0.9
Alaska-Inkutzk	CLB	IKT-04-3	duplicate	0.29	Not available	25.364	0.92
Alaska-Inkutzk	CLB	IKT-04-4		0.5	86.041	232.207	2.58
		missing					
Alaska-Inkutzk	CLB	IKT-04-5		0.15	32.872	21.066	4.24
Alaska-Inkutzk	CLB	IKT-04-5		0.17	Not available	1.775	0.78
Alaska-Inkutzk	CLB	IKT-04-6		0.85			
Alaska-Inkutzk	CLB	IKT-04-6		0.32			
Alaska-Inkutzk	CLB	IKT-04-7		0.17			
Alaska-Inkutzk	CLB	IKT-04-7		0.28			
Alaska-Inkutzk	CLB	IKT-04-8		0.74			
Alaska-Inkutzk	CLB	IKT-04-8		0.46			
Alaska-Inkutzk	CLB	IKT-04-9		5.248	3.178		
Alaska-Inkutzk	CLB	IKT-04-9	duplicate	0.67			
Alaska-Inkutzk	CLB	IKT-04-9		0.35			
Alaska-Inkutzk	CLB	IKT-04-9	duplicate	0.35			
Alaska-Inkutzk	CLB	IKT-04-10		0.76			
Alaska-Inkutzk	CLB	IKT-04-10		0.27			






# Finding the Differentiator

*The Flux:*

*Laxman: Sodium Sulfate  
(Na<sub>2</sub>SO<sub>4</sub>)*

*Common Practice:  
potash (K<sub>2</sub>CO<sub>3</sub>)*

*Laxman's glass would  
have a lower K to Ca  
ratio than potash-based  
glasses.*




**Thermo**  
S C I E N T I F I C




## *Thermo Scientific NITON Analyzer*

- Helium (He) purge feature of these analyzers allows for analysis of light elements Mg, Si, P, S, K, Ca.*
- Nondestructive Analysis*
- Short analysis time*
- Spectra and (semi)quantitative data*




*Figure 1.* Light energy spectra collected from glass samples from Castle Hill, Alaska. Potassium and calcium are labeled, showing an average K to Ca ratio between 1:2 and 1:1.







*Figure 2.* Light energy spectra collected from glass samples from Yakutat, New Russia. Potassium and calcium are labeled, showing an average K to Ca ratio between 1:1 and 2:1.




*Figure 3.* Light energy spectra collected from glass samples from Irkutsk, the site of Erik Laxman's glass factory. Potassium and calcium are labeled, showing an average K to Ca ratio between 1:2 and 1:1, with one unusual outlier with a much greater concentration of K.




*Figure 4. Light energy spectra from both Irkutsk (green) and Yakutat (blue) overlaid.*



**Figure 5.** Light energy spectra from both Castle Hill (red) and Yakutat (blue) overlaid.




*Figure 6. Light energy spectra from both Castle Hill (red) and Irkutsk (green) overlaid.*



**Figure 8.** Light energy spectra collected from glass samples from Castle Hill. Sulfur peak labeled.

**Figure 7.** Light energy spectra collected from glass samples from Irkutsk. Sulfur peak labeled.



# Discussion and Conclusions

- The Irkutsk glass samples have a lower K to Ca ratio as compared to the Yakutat samples, and a similar ratio as compared to the Castle Hill samples. We will need to look more closely at quantitative data, rather than spectra to better establish these ratios.
- The Irkutsk samples have a significantly lower concentration of potassium than both the Yakutat and Castle Hill samples. Further study is necessary to determine the degree of significance of ratio over concentration, or vice versa.
- Overall, the results lead us to believe that with better method and a more statistically significant number of samples, the handheld XRF is a viable tool for this research.
- Several follow-up steps need to be taken for more conclusive results:
  - A more statistically significant number of samples from each location needs to be analyzed
  - Samples need to be better matched (color, use, etc.) for comparison
  - Comparison of a greater number of elements needs to be undertaken
  - Ratios need to be more precisely established by using the available (semi)quantitative data.

# References

Bychkov, Oleg V.

- 1997 The Origin of Colonial Glass Production in Irkutsk: Research Perspectives. In *The Archaeology of Russian Colonialism in the North and Tropical Pacific*, edited by Peter R. Mills and Antoinette Martinez. *Kroeber Anthropological Society Papers* 81:42-49.

Crowell, Aron L.

- 1997 *Archaeology and the Capitalist World System: A Study from Russian America*. Plenum Press, New York and London.

McMahan, J. David

- 2004 *Archaeological Investigations at Castle Hill: the Early 19th Century Colonial Capitol of Russian-America, Sitka, Alaska*. Manuscript for publication submitted to Dr. Alexander Artemiev, Institute of History, Archaeology and Ethnography of the Far Eastern Branch of the Russian Academy of Sciences.

- 1999 Archaeology at Alaska's Castle Hill, the Colonial Capital of Russian America. In *Unseen Treasures: Imperial Russia and the New World*, pp. 17-22. American-Russian Cultural Cooperation Foundation, Washington, D.C., and the Russian State Historical Museum, Moscow.